

AZ 可以提供什么样的 PH 探针?

1. AZ 制造的所有 pH 探头均为玻璃球电极、单参比接界，非再填充型，参考系统为 Ag/AgCl，参考电解质为 KCl。

2. AZ 目前有 9 个型号的 pH 探头，如下：

型号照片	规格	可搭配之仪表型号
 <p>VZ86P2AZ</p>	<p>pH 范围：0.00 - 14.00 测量温度：0 至 80°C 温度补偿：不适用 直径：12 毫米 轴长：140 毫米 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>
 <p>VZ86P3AZ</p>	<p>pH 范围：0.00 - 14.00 工作温度：0 至 80°C 温度补偿：是 （热敏电阻） 直径：12 毫米 轴长：140 毫米 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>
 <p>VZ86P4AZ</p>	<p>pH 范围：0.00 - 14.00 专为低离子强度而设计！ 工作温度：0 至 80°C 温度补偿：不适用 直径：12 毫米 轴长：140 毫米 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>
<p>玻璃</p>  <p>VZ86P6AZ</p>	<p>pH 范围：0.00 - 14.00 专为低离子强度而设计！ 全玻璃球形 pH 探头 工作温度：0 至 80°C 直径：12 毫米 轴长：145 毫米 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>

 <p>VZ86P7AZ 尖头</p>	<p>pH 范围：0.00 - 14.00 专为低离子强度而设计！ Sharp Glass Tip pH 探头 工作温度：0 至 80° C 直径：12 毫米 轴长：143 毫米 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>
 <p>VZ86P1AZ 平面</p>	<p>pH 范围：2.00 - 12.00 专为低离子强度而设计！ 平板计算机表面 pH 探头 工作温度：0 至 80° C 直径：12 毫米 轴长：134 毫米 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>
 <p>VZ86P11AZ</p>	<p>pH 范围：0.00 - 14.00 专为低离子强度而设计！ 6 mm 玻璃球形 pH 探头 工作温度：0 至 80° C 直径：6 毫米 轴长：待定 连接器：BNC 型 电缆长度：1238mm</p>	<p>86501,86502,86504,86505, 86551,86552,86554,86555, 9861,9661,8651,8601</p>
 <p>VZ863PAZ</p>	<p>pH 范围：2.00 - 12.00 专为低离子强度而设计！ 工作温度：0 至 60° C 温度补偿：是 （热敏电阻） 直径：22 毫米 轴长：150 毫米 连接器：IP67 针型 电缆长度：1000 毫米</p>	<p>84051, 86031, 86021</p>

 <p>VZ861PAZ</p>	<p>pH 范围: 2.00 – 12.00 专为低离子强度而设计! 工作温度: 0 至 60°C 温度补偿: 是 (热敏电阻) 直径: 22 mm 轴长: 150 mm 连接器: IP67 针型 电缆长度: 1000 mm</p>	<p>8405, 8603, 8602</p>
---	--	-------------------------

对于台式计、便携式 pH 计，AZ 可以提供以下准确的 pH 探头，

1. 86P3: ATC pH 探头
2. 86P2: 不带 ATC 功能的 pH 探头
3. 86P4: 低离子强度的 pH 探头，无 ATC 功能
4. 86P6: 全玻璃 pH 探头，不带 ATC 功能
5. 86P7: 无 ATC 功能的尖头 pH 探头
6. 86P1: 不带 ATC 功能的平面 pH 探头
7. 86P11: 6mm 玻璃 pH 探头，不带 ATC 功能

对于 AZ Combo 水质测试仪(例如: 86031、86021、84051), 其 pH 探头为“863P”。

For the Benchtop Meter, Portable pH Meter, AZ could provide the following accurate pH probe,

1. **86P3**: ATC pH probe
2. **86P2**: pH probe without ATC function
3. **86P4**: pH probe for Low Ion Strength without ATC function
4. **86P6**: Full glass pH probe without ATC function
5. **86P7**: Sharp tip pH probe without ATC function
6. **86P1**: Flat surface pH probe without ATC function

7. **86P11**: 6mm glass pH probe without ATC function

For the AZ Combo Water Testing Meters (ex: 86031, 86021, 84051), their pH probe is "**863P**".